
Liquid handling

LIQUID HANDLING

LIQUID HANDLING

LIQUID HANDLING

**DOSIPET**

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DOSIPET

Dosipet is a series of adjustable volume repetitive dosing devices. They are made of high quality materials ensuring considerable mechanical and chemical resistance. They are available in three different versions, depending on the reagents to be dispensed. Possibility to choose between:

DOSIPET B

For universal use; ideal for use with reagent of average aggressiveness.

DOSIPET S

Useful for dosing concentrated acids, polar and non-polar solvents.

DOSIPET HF

Especially developed for hydrofluoric acid dosing. They are supplied complete with tubes for filling/discharging, as well as with a set of adaptors: for Dosipet 2.5/5 and 10 ml adaptors A25, A28, A45, S40; for Dosipet 25/50 and 100 ml adaptors A28, A32/45, A40.

Technical specifications common to the three series:

- precise and quick volume adjustment using the analogue slider
- precise graduated scale to ensure reproducibility of the dispensed volumes
- no loss of reagent when eliminating air or filling
- dispensing without bubbles forming
- the piston completely transfers the liquid with every movement
- total recovery of the liquid in the discharge tube by turning the delivery tip 180°
- autoclavable at +121°C for 20 minutes
- EC conformity

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ALKALI	INORGANIC SOLUTIONS	ACIDS	ORGANIC SOLVENTS		HYDROFLUORIC ACID (HF)
			POLAR	NON-POLAR	
Dosipet					
			Dosipet S		
			Dosipet HF		

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DOSIPET B

Series that can be used with every water solution type (except those containing hydrofluoric acid), marked by a high degree of convenience and reliability. It is equipped with a 3.3 borosilicate glass cylinder and PTFE piston

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DOSIPET S

High quality components and assembly make this series the ideal choice for measuring out organic polar and non-polar solvents (see technical manual), including halogenated solvents, and concentrated mineral acids, including particularly volatile acids such as nitric and trifluoroacetic acids. In addition, they can be used with peroxides, hydrochloric acid and sulphuric acid. It comes with a borosilicate glass piston. The reagents come into contact only with parts made of PTFE.

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DOSIPET HF

These dosing devices are designed for specific use with hydrofluoric acid solutions. The internal components that come into contact with these solutions are made with high quality and high resistance materials (PTFE and platinum-iridium), completely inert with hydrofluoric acid. Complete with PTFE piston.

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DOSILAB

Dosilab is a new series of adjustable volume bottle top dispenser for dosing chemicals with precision, accuracy and inexpensiveness. It is totally safe for operators and environment. The six models of the new line are simple to use and easy to disassemble to perform common cleaning and maintenance operations; autoclavable without need for disassembly. Volume adjustment is made precise and quick by using the analogue slider. A careful choice of raw materials, as well as several quality controls carried out during production process guarantee Dosilab's high precision and accuracy. Every single appliance is individually calibrated according to the ISO 8655 standards and it is supplied with its own calibration certificate. The supply includes five adaptors for use with the most common reagent bottles (28, 32, 36, 40 and 45 mm) and a practical device for easy recalibrations that can be directly performed in laboratory, in conformity with the requirements of the GLO/ISO standards.

TECHNICAL SPECIFICATIONS

- Completely autoclavable at +121°C without need for disassembly.
- PTFE piston with silicone O-ring ensuring a high degree of accuracy.
- Internal cylinder in borosilicate glass protected by PP sheath.
- Valve housing: PTFE body in a single piece ensuring high chemical resistance.
- Accuracy & Reproducibility: according to the ISO 8655 standards
- Compatibility: excellent with all reagents including organic solvents and concentrated acids, except hydrofluoric acid (HF).
- Recalibration: simple and quick with its special accessory.

Comparison of the specifications required by the ISO 8655 standards and the DOSILAB specifications

The following table certifies that the DOSILAB specifications fully meet the requisites required by ISO 8655, international standard for piston-driven volumetric appliances.

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All the internal components that during use are in contact with liquid solutions are made of PTFE, FEP and borosilicate glass.

They therefore possess excellent compatibility with all reagents (even organic solvents and concentrated acids), except for hydrofluoric acid.

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Alkali	Inorganic Solutions	Acids	Organic Solvents		Hydrofluoric Acid (HF)
			Polar	Non-polar	
Dosilab					

Art.	Mod.	Cap. ml	Subdiv. ml	Prec. ml	CV ml	Filettatura GL / Thread GL / Filetage GL
83082	Dosilab 2	0,2 - 2,5	0,05	0,015	0,005	30
83083	Dosilab 5	0,5 - 5	0,10	0,025	0,005	30
83084	Dosilab 10	1 - 10	0,20	0,050	0,010	30
83085	Dosilab 30	2,5 - 30	0,50	0,150	0,030	30
83086	Dosilab 60	5 - 60	1,00	0,300	0,060	30
83087	Dosilab 100	10 - 100	2,00	0,500	0,100	30

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TECHNOTRATE DIGITAL BURETTE

Technorate is a digital burette ideal for exact titrations directly from bottles, both in the laboratory and in the field. Developed in compliance with EN ISO 9001, it is EC certified. It has a plastic body and a large liquid crystal display also for continuous display of the dispensed volume. It adequately allows dosage volumes of up to 50 ml with 0.01 ml increases to be covered and the dosage velocity can be selected from fast mode to drop-by-drop mode. Technorate can store up to 1,000 dosage and titration values entered by the user, which can be downloaded onto computers or printers using the specially provided interface.

Simple calibration

Permanent calibration refers to water at +20°C. In any case, when using liquids with a different viscosity/density, Technorate can easily and quickly be recalibrated by pressing the dedicated button.

Practical and quick preselection of the titration volumes

Titration volume preselection is made particularly simple by pressing only the dedicated "DOS" button, looking at the display, and with just a few turns of the handle the burette cylinder is filled.

The special memory button ("MEM") helps you avoid complex calculations.

Easy and safe

By turning 180°, the delivery tip makes all the liquid in the discharge tube flows inside the bottle and prevents the spilling of inappropriate liquids. The components that come into contact with the reagents are made of high resistance plastic.

Technorate is tested one by one at the end of production, has quality certification and an individual serial number.

It is supplied in series with A25, A28, A32/45 and S45 threaded adaptators and PTFE aspiration, discharge and recirculation tubes

TECHNICAL SPECIFICATIONS:

- auto-off function after 3 minutes of non-use
- use temperature from +10°C to +40°C
- transport and storage no warmer than +60°C
- rechargeable battery: the intermittent indicator inside the display informs of the need for recharge without jeopardising the titration in progress
- precision mechanism for simple filling and elimination of air without reagent loss problems
- dispensing without bubbles forming
- battery charger: 100-240V/ 50-60Hz

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VARIABLE VOLUME AUTOCLAVABILE MICROPIPETTES - PLURIPET

Eight variable volume single-channel micropipettes make up the PLURIPET line equipped with an accurate and precise pipetting mechanism. They cover the range of capacities from 0.2 µl to 10,000 µl. The PL2 and PL10 models are ideal for molecular biology techniques. The PL20, PL100, PL200 and PL1.000 models are useful for dispensing small volumes of aqueous fluids of moderate density and viscosity. The PL5.000 and PL10.000 macropipettes are ideal for dispensing large volumes of rather thick and viscous aqueous fluids. They are equipped with a special filter that protects the internal piston by increasing its lifetime and preventing contaminations. They are supplied complete with a recalibration device and adhesive support for attachment to the laboratory bench. Individually calibrated with their certificates included in the package; EC/IVD-compliant.

TECHNICAL SPECIFICATIONS

- the handle connection ensures a more comfortable and non-fatiguing grip
- the volume can be selected by turning the emptying button or turning the cylinder supplied for the purpose
- the large display ensures simple reading of the set volume value
- the button of each micropipetta has a specific colour so identification of the usable volume interval is simple
- the tip ejector is adjustable in height to house a large variety of tips and is can be disassembled to simplify cleaning operations
- the recalibration device ensures quick and precise recalibrations
- the robust structure of the PVDF body ensures high mechanical and chemical resistances
- totally autoclavable at +121°C for 20 minutes
- the materials used for the structure are resistant and durable, and permit exposure to U.V. rays
- in accordance with EN IS 8655

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FIXED VOLUME AUTOCLAVABLE MICROPIPETTES - FIXOPET

The FIXOPET line is formed by eleven single channel micropipettes with fixed volume of high precision. The capacities available range from 5 µl to 1,000 µl. Each of them is equipped with adjustable tip ejector in order to eliminate any risk of contamination. The structural features of these micropipettes considerably reduce the force necessary for pipetting. Ideal for clinical and educational uses, compliant with EC/IVD.

TECHNICAL SPECIFICATIONS

- The special handle connection ensures a more comfortable and non-fatiguing grip
- the rotary device for the tip ejection does not create strain on the tendons of the thumb
- the removable ejector allows it to be used with test tubes having a particularly narrow opening
- the integrated recalibration system ensures quick and precise recalibrations
- the robust structure of the body made of PVDF ensures high mechanical and chemical resistances
- the materials used for the structure are resistant and durable, and permit exposure to U.V. rays
- totally autoclavable at +121°C for 20 minutes
- in accordance with EN IS 8655

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SPARE PARTS SET AND MICROCHANNEL MICROPIPETTES STAND

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Bench holder in transparent plastic material for 3 up to 6 Kartell micropipettes (Pluripet, Pluripet II, Fixopet, Labofix) or other brands. Strong structure, easy cleaning, wide base allowing great stability. Radiation sterilizable.

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LIQUID HANDLING

MULTICHANNEL VARIABLE VOLUME MICROPIPETTES - MULTIPET

Multichannel micropipettes with air cushion ideal for dispensing in multiwell plates for ELISA dosages, molecular biology and DNA amplifications. The special PP/PVDF structure makes them lightweight and comfortable. The pipetting mechanism makes the work very precise and the adjustment, which can be carried out with just one hand, is simple and progressive. The exclusive operation system reduces the force necessary for executing the pipetting and suction; it considerably reduces the force necessary for inserting and ejecting tips.

TECHNICAL SPECIFICATIONS

- The volume can be easily selected by turning the emptying button or turning the cylinder supplied for the purpose
- the large display ensures easy reading of the set volume value
- the MULTIPET cones are independent, so the required force for securing the 8 or 12 tips is uniformly distributed: minimum pressure suffices for inserting and ejecting them. In addition, the movement of the cones allows the micropipette to oscillate during the loading operations
- the lower part can freely turn 360°, the ejector is easily adjustable to house a large variety of tips and can be removed for cleaning
- they are totally autoclavable at +121°C for 20 minutes, without having to recalibrate them
- the materials used for the structure are resistant and durable, and permit exposure to U.V. rays
- individually calibrated with their certificates included in the package; EC/IVD-compliant
- **in accordance with EN IS 8655**

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FIXED VOLUME AUTOCLAVABLE MICROPIPETTES - LABOFIX

This line of micropipettes is designed to meet the increased demand for fixed volume micropipettes that combine a good quality level with a certainly interesting price. Available ranging from the volume of 5 µl up to the volume of 1,000 µl.

Ideal for applications in Clinical Diagnostics, Molecular Biology, Cellular Biology and Analytical Chemistry and in all those fields in which a fixed dosage carried out with good precision and accuracy is necessary.

The particularly lightweight body makes the handle comfortable. Individually calibrated according to the ISO8655 standards, the new international standards for piston-driven volumetric appliances; EC/IVD-compliant.

TECHNICAL SPECIFICATIONS

- made of high quality plastic
- equipped with tip ejector and stainless steel piston for use prolonged over time
- the handle connection ensures a more comfortable and non-fatiguing grip
- the tip holder cone can house a wide variety of tips; it has a fine shape to make getting inside test tubes or containers having particularly narrow openings or lengthened shapes easier
- the robust structure of the body ensures high mechanical and chemical resistances, making them totally autoclavable at +121°C for 20 minutes

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VARIABLE VOLUME AUTOCLAVABLE MICROPIPETTES - PLURIPET II

Line made up of eight models having air displacement and equipped with coaxial ejector. Useful for applications in Molecular Biology, Microbiology, Immunology, Analytical Chemistry, etc. They ensure a high degree of precision and accuracy combined with a lower cost. The particularly light and balanced body makes the handle comfortable. They cover the range of capacities from 0.2 µl to 10,000 µl. The Pluripet II 5,000 and Pluripet II 10,000 macropipettes are ideal for dispensing large volumes of rather thick and viscous aqueous fluids. They are supplied complete with adhesive support for attachment to the laboratory bench. Individually calibrated according to the ISO8655 standards, the new international standards for piston-driven volumetric appliances; EC/IVD-compliant.

TECHNICAL SPECIFICATIONS

- made of high quality plastic and equipped with steel piston for prolonged use over time
- the handle connection ensures a more comfortable and non-fatiguing grip
- the desired volume can be set by turning the emptying button; it proves easy and safe thanks to the snap selection
- the large display ensures simple reading of the set volume value
- the button of each micropipette has a specific colour so identification of the usable volume interval is simple
- the tip holder cone can house a wide variety of tips
- the recalibration device ensures quick and precise recalibrations
- the robust structure of the body ensures high mechanical and chemical resistances, making them totally autoclavable at +121°C for 20 minutes

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MINIFIX

An innovative product line including ten fixed-volume air displacement mini-micropipettes. Fully autoclavable, they are designed for fast, repetitive and precise operation and are competitively priced. Very compact in terms of size and weight, each mini-micropipette is 13 cm long and weighs approximately 10 grams. Designed for routine or biological hazard work, in situ analyses and sampling, diagnostic kits and instruction. Kartell tips fit on the conical shaft as well as most other commercially available tips. The plunger button has a different color based on dosing capacity. Individually calibrated in compliance with ISO8655, the new international standard for piston-operated volumetric apparatus. Carries the CE conformity mark.

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LIQUID HANDLING

SEROPET

Seropet was developed for aspirating liquids by glass or plastic pipettes ranging from the volume of 0.5 ml up to 100 ml. The particularly light structure ensures comfortable use even when used for a long time and its particular configuration also makes it ideal for underhood uses. The buttons are very easy to access for choosing different modes depending on the volume pipette being used and viscosity of the liquid to aspirate. It is protected from crossed contamination in the event of excessive aspiration by the filter on the aspiration device and the safety valve. It is supplied complete with battery charger, protective PTFE hydrophobic filter (0.2 µm) and an adhesive support for attachment to the laboratory bench. EC/IVD-compliant.

TECHNICAL SPECIFICATIONS

- two different work modes: quick and normal, with continuous or pulsating adjustment
- the pipette housing is autoclavable
- supplied with safety valve to block aspiration excesses
- equipped with autoclavable protective hydrophobic anti-aerosol filter
- the special parking fins make it possible to set it down on the bench to prevent contamination of the pipette being used
- highly sensitive internal valves to get precision also with small volumes
- the NiMH battery supplied has an autonomy of up to 8 hours with continuous work and up to 1,000 charge cycles
- the indicator light warns when the battery charge is about to become exhausted

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**LABOFILL - CORDLESS
ELECTRONIC PIPETTE FILLING
DEVICE LABOFILL**

A new pipette filling device was developed to complete Kartell entry level line of liquid handling products. This device is used for filling and dispensing liquid in glass or plastic pipettes, both graduated and volumetric, up to 100 ml. LaboFill is a safe, efficient and economical instrument for liquid filling and dispensing. Due to its design is also comfortable and easy to use.

LaboFill is CE marked, is individually tested after production process and is calibrated according to ISO 8655 standards.

PRODUCT FEATURES:

- Red and green L.E.D.s: the red light indicates low battery charge; green L.E.D. indicates charging status.
- Thumb wheel is provided at the back of the instrument to increase or reduce pump speed according to user's needs.
- High aspirating and flow-out speed: 25 ml in only 5 seconds.
- Silicon universal pipette housing can accommodate standard pipettes from 0.1ml to 100 ml. Inside collet knurlings facilitate tight fit even with wet pipettes.
- Comfortable handgrip ensures optimal grip and prevent user from tiredness
- 0.2µm PTFE hydrophobic filter is provided in the collet assembly to avoid aerosol formation and liquid entering the instrument body.
- Rechargeable battery Ni-MH 3.6 provided with the instrument, ensures cordless work for up to eight hours continuing after recharging.
- Completely battery recharging in 5 hours.
- Adaptor compatibility with 220 V and 110 V

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